

She's on Top of the World

Andraea LaVant has places to go and people to see! With safe seat elevation on her 4Front® 2, Andraea is easily visible on the streets of Tempe, Arizona.

 [CLICK HERE](#) to watch the video.



Center This

The basic center mount foot platform is now available!


FEATURES:

- 65 to 100 degrees of angle adjustment
- Height adjustment: 8-10", 10-14", 14-19"
- Angle adjustable foot plates
- Mounts to inside of side rail (TB Flex, Q captain's seat) or depth-adjustable bracket (TB3 seating, iLevel)

Save Time

Here are four ways the PQ-365 app saves you time:

- Serial number/order tracking quick search: find individual components, no more scrolling
- Fault code search: look up codes while in the field
- Brochures/spec sheets: Check product specs while on the go
- Installation instructions: access step-by-step guides for parts, products and accessories



SCAN NOW TO INSTALL

PQ-365 is the one-stop app for wherever you go, whenever you need it.



Quantum® news updates are designed to keep NSM ATPs current on everything new. If you have any questions, please call Quantum Quick Connect at **877-800-7808**.