




# Q-LOGIC 3<sup>e</sup>



## SEPTEMBER UPDATE

News you can use in 30 seconds or less!



### It's Only Logical

Updates to Q-Logic 3e are now available.

All-new features:

- Industry-first, magnetic charging port with adapter can be used with any current XLR charger
- Magnetic wall holder keeps charger plug within reach and off the floor
- Integrated puddle light provides greater visibility
- Q-Logic 3e comes standard on the 4Front® 2 and 4Front® 2 HD

[CLICK HERE](#) to watch the video.

### Make Your Voice Heard

The Centers for Medicare & Medicaid Services has opened a 30-day public comment period on national coverage analysis for power seat elevation systems, such as iLevel®, as an accessory to power wheelchairs (Group 3).

[CLICK HERE](#) to leave a comment.



### Consider This

When choosing a power chair for your client, ask yourself these questions:

- Can they perform independent pressure relief?
- Can they complete MRADLs safely, independently and in a timely manner?
- Do they require alternative drive controls?
- Does the home have adequate access for maneuvering a chair?



[CLICK HERE](#) to read more.

### On the Road

Need an answer on the road? The PQ-365 app is a great solution. Look up product specs, IPBs, or serial numbers.

**SCAN NOW TO INSTALL**

PQ-365 is the one-stop app for wherever you go, whenever you need it.



FDA Class II Medical Device

### A Ride Like No Other

The Jazzy® EVO 613 Series is narrow, maneuverable and delivers a smooth, comfortable ride.

- Active-Trac Suspension combines with Mid-Wheel® 6 Technology
- Active front casters lift when climbing a threshold or obstacle
- Casters also drop down for added stability when descending

[CLICK HERE](#) to read more.