



TRU BALANCE[®] 4

POWER POSITIONING SYSTEMS

Now Available on **Stretto, Edge[®] 3** and **4Front[®] 2**

30°
20°
10°



Power to the Position

Now available for ordering!
TRU-Balance[®] 4 Power Positioning Systems with Anterior Tilt.

- Anterior tilt options: 10°, 20° and 30°
- Up to 8 customizable memory seating positions
- Articulating Foot Platform lowers to assist with transfers
- Available on the Stretto, Edge[®] 3 and 4Front[®] 2

AFP articulation range: 0-7.25"

AFP elevation range: 85-165°



 [CLICK HERE](#) to watch the video.



We Got the Need for Speed

Our ambassadors tell us how fast they can go while elevated at iLevel[®].

 [CLICK HERE](#) to watch the video.




Charge!

Proper charging habits are key to ensure new or experienced power chair users get the most out of their battery life.

- Chair used every day - charge 8-14 hours every night
- Chair used less frequently - charge 8-14 hours 2x week
- Do not charge longer than 14 hours

 [CLICK HERE](#) to learn more.






Break it Down

Illustrated Parts Breakdowns (IPBs) simplify the parts search for Pride[®] and Quantum[®] products. Download the PQ-365 app today.

SCAN NOW TO INSTALL

PQ-365 is the one-stop app for wherever you go, whenever you need it.



Quantum[®] news updates are designed to keep NSM ATPs current on everything new. If you have any questions, please call Quantum Quick Connect at **877-800-7808**.